

Test Sequence

IBP Monitoring – IEC 60601-2-34 (v1.0)

Template	IBP Monitoring – IEC 60601-2-34
Version	1.0
Created	2025-10-02 01:08
Created By	Seed

Test Steps

#	Step / Instruction	Type	Limits	Unit	Actual / Pass
1	Arterial Pressure Accuracy (120/80) Simulate arterial pressure waveform; record systolic/ diastolic. Ref: IEC 60601-2-34 §201.12.4.102	Measurement	120 mmHg ±4	mmHg	-----