

Test Sequence

Patient Hoist – EN ISO 10535 (v1.0)

Template	Patient Hoist – EN ISO 10535
Version	1.0
Created	2025-10-02 01:08
Created By	Seed

Test Steps

#	Step / Instruction	Type	Limits	Unit	Actual / Pass
1	Safe Working Load (Label Check) Confirm SWL label and verify test load setup. Ref: EN ISO 10535 §4.2.2	Measurement	200 kg ±0	kg	_____
2	Emergency Lowering Function Test emergency lowering under load simulation. Ref: EN ISO 10535 §4.3.4	Yes/No Check	-	-	_____