

Test Sequence

Ambulance Trolley/Stretchers – EN 1865 (v1.0)

Template	Ambulance Trolley/Stretchers – EN 1865
Version	1.0
Created	2025-10-02 01:08
Created By	Seed

Test Steps

#	Step / Instruction	Type	Limits	Unit	Actual / Pass
1	Locking & Height Adjustment Check safe locking at all heights; smooth operation. Ref: EN 1865-1	Yes/No Check	-	-	-----
2	Load & Stability Check Verify safe working load and tip resistance per site policy. Ref: EN 1865-1 / EN 1789	Yes/No Check	-	-	-----