

Test Sequence

NIBP Monitoring – IEC/ISO 80601-2-30 (v1.0)

Template	NIBP Monitoring – IEC/ISO 80601-2-30
Version	1.0
Created	2025-10-02 01:07
Created By	Seed

Test Steps

#	Step / Instruction	Type	Limits	Unit	Actual / Pass
1	BP Accuracy 120/80 mmHg Simulate 120/80 via NIBP tester; record reading. Ref: ISO 80601-2-30 §201.12.1	Measurement	120 mmHg ±5	mmHg	_____
2	Diastolic Accuracy 80 mmHg Verify diastolic reading. Ref: ISO 80601-2-30 §201.12.1	Measurement	80 mmHg ±3	mmHg	_____